

Art-Can Safety

Safety Arterial Cannula



Art-Can Safety

Safety Arterial Cannula



On / OFF
Switch



Needle Safety
Guard

Advantages: Safety and Comfort

- **Specially designed ON-OFF switch**
Prevents backflow and exposure of blood to healthcare professional and provides safety from air embolism and contamination
- **Needle safety guard**
Automatically covers the needle's sharp bevel after withdrawal of needle from the hub minimizing the risk of Needle Stick Injury
- **Bio-Compatibility**
Advance Polyurethane catheter material ensures better bio-compatibility and longer indwelling time
- Also available in FEP catheter material
- Also available in standard version (without needle safety guard)
- Sizes: 20G & 22G

Applications: Suitable for arterial pressure monitoring and serial arterial blood gas analysis

Sales & Marketing Office :

Poly Medicure Limited
232-B, 3rd Floor, Okhla Industrial Estate,
Phase-III, New Delhi-110020, India
Tel:+91-11-33550700, 47317000
Fax:+91-11-26321894/1839
Customer Care No. : +91-129-4287053
info@polymedicure.com

Europe Representative Office :

Via Pietro. Verri 1
Business Center BOMA,
Build 1C 46100, Mantova. Italy
polymed.eu@polymedicure.com

EC REP

Obelis s.a
Boulevard Général Wahis 53
1030 Brussels, BELGIUM
Tel. : + (32) 27325954
E-mail : mail@obelis.net